


Analog Sensor Extension Cable

(Order Code EXT-BTA)



This cable can be used with any Vernier sensor that connects to an analog channel of LabQuest[®] 2, LabQuest, LabPro[®], TI-Nspire[™] Lab Cradle, CBL 2[™], or CBL[™] such as a pH Sensor or Dual-Range Force Sensor. Connect the plug end of the cable to the appropriate port on the lab interface and the socket end to the sensor. If additional length is needed, you can connect additional extension cables. Please contact us if you have any other questions about using this cable.



Measure. Analyze. Learn.[™]

Vernier Software & Technology

13979 S.W. Millikan Way • Beaverton, OR 97005-2886

Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440

info@vernier.com • www.vernier.com