

# Digital Sensor Extension Cable

(Order Code EXT-BTD)



This cable can be used with any Vernier sensor that connects to a digital channel of LabQuest<sup>®</sup> 2, LabQuest, LabPro<sup>®</sup>, TI-Nspire<sup>™</sup> Lab Cradle, CBL 2<sup>™</sup>, or CBL<sup>™</sup> such as Motion Detectors and photogates. Connect the plug end of the cable to the appropriate port on the lab interface and the socket end to the sensor. If additional length is needed, you can connect additional extension cables. Please contact us if you have any other questions about using this cable.



Measure. Analyze. Learn.™

**Vernier Software & Technology**

13979 S.W. Millikan Way • Beaverton, OR 97005-2886

Toll Free (888) 837-6437 • (503) 277-2299 • FAX (503) 277-2440

info@vernier.com • www.vernier.com